# CEREC Ortho SW 2.0.1

## Update Description (EN)

### Installation notes:

- Perform a backup of all patient data before starting the installation process.
- Do not perform an installation on several acquisition units or computers at the same time if these are within the same network.
- Before starting the installation process, please check the patient database directory in the software (Configuration/Settings/Patient Database). The path should always point to a directory, never to the root of a drive, i. e. never to "D:\", but rather to "D:\Data".

### Hardware support note:

CEREC Ortho SW 2.0.1 supports Omnicam (with guided scanning) and Primescan (no guided scan required) scanners, but requires the newer AC models with touchscreen and Windows 10 to run. An upgrade version for existing Omnicam installations will follow at a later date.

### Changes made to previous version 2.0:

- Application Setup
  - Fixed a bug which prevented the version from being installed on acquisition centers with the newest hardware configuration.
- ACQUISITION Phase:
  - Fixed memory leaks
  - Improved memory handling after deleting image catalogs
  - Progressive image reduction to improve performance during longer scans
  - Shortened delay when resuming scan process
  - Significantly reduced application crash probability
- Model Analysis:
  - Fixed L/R labels to patient orientation
  - Bugfix for missing age for some languages in PDF export
  - Added missing Moyers and Nance values to PDF export

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