Sidexis 4

- Installationsvoraussetzungen
- Installation Requirements
- Conditions d'installation
- Requisitos de instalación
- Requisiti di installazione

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This document covers the most important conditions that must be met for the installation of Sidexis 4 version 4.3.1 in practice.

The following aspects need to be clarified beforehand with the customer:

- PC hardware
- Devices put to practical use
- Sidexis updates

PC hardware

PC server

Requirements for installation of Sidexis 4:

It is assumed for each operating system to be of the current patch state. There must be an internet connection running during installation.

Sidexis 4 requires an MS SQL Server. Sidexis 4 is compatible with the MS SQL Server versions 2014 and 2017. The MS SQL Server 2014 Express is brought along and installed by the Sidexis 4 setup, if a compatible version is not already installed.

Security note: Please keep your operating system and other Microsoft components up-to-date with security updates even after installation.

	Minimum requirements	Recommended requirements
Operating system	Windows 7 Pro SP1 (64 bit) ¹ Windows 8.1 Pro (64 bit) ¹ Windows 10 Pro (64 bit) ^{1, 2} Windows Server 2008 R2 Windows Server 2012 R2 Windows Server 2016 Windows Server 2019	Windows Server 2008 R2 Windows Server 2012 R2 Windows Server 2016 Windows Server 2019
RAM	≥ 4 GB	≥ 8 GB
CPU	≥ 2 GHz DualCore	≥ 2.3 GHz QuadCore processor with 64bit (x64)
Hard disk	> 675 GB	> 1 TB
	Sufficient disk space must be ensured even during operation	
Graphic card	≥ 512 MB	≥ 1GB
	Graphic card with dedicated graphic memory or Intel Onboard Graphic ³	
DirectX	DirectX 9.0c	DirectX 10
	with WDDM driver 1.0 or higher	

1. Using these systems for server operation is not recommended, as they are designed and scaled up by Microsoft for this purpose.
2. Tested: Version 1903 Build 18362.175.
3. The Intel Onboard Graphic has no graphic memory of its own and uses the PC RAM.

In the case of installation on a domain controller, please observe the notes regarding this in the Sidexis 4 Service Manual.

Workstation PC

Requirements for installation of Sidexis 4: It is assumed for each operating system to be of the current patch state. There must be an internet connection running during installation.

	Minimum for 2D-stations	Minimum for 3D-stations
Operating system¹	Windows 7 Pro SP1 (32 or 64 bit) Windows 8.1 Pro (64 bit) Windows 10 Pro (64 bit) ²	Windows 7 Pro SP1 (64 bit) Windows 8.1 Pro (64 bit) Windows 10 Pro (64 bit) ²
RAM	≥ 4 GB	≥ 8 GB
CPU	≥ 2 GHz DualCore	> 2.3 GHz QuadCore processor with 64 bit (x64)
Hard disk	≥ 5 GB	≥ 5 GB
Graphic card	≥ 512 MB Graphic card with dedicated graphic memory or Intel Onboard Graphic ³	≥ 1GB
DirectX	DirectX 9.0c	DirectX 10 with WDDM 1.0 or higher driver

1. For 64-bit operating systems, installation under Bootcamp is approved
2. Tested: Version 1903 Build 18362.175.
3. The Intel Onboard Graphic has no graphic memory of its own and uses the PC RAM.

In order for the interaction with the 3D rendered volume to be reliably "judder-free", graphics cards are recommended with at least the following Passmark GPU benchmark values:

NVIDIA: Passmark > 1000

AMD: Passmark > 1500

Intel Onboard: Passmark > 540

The Passmark values can be checked, if the name of the graphics card is known, at www.videocardbenchmark.net.

Screen

	Minimum requirements	Recommended requirement
Resolution	1280 x 1024 pixels	1600 x 1200 pixels

Sidexis 4 practical devices

All devices in the current imaging systems portfolio are supported

supported	no longer supported
Xios	Sidexis IO sensors
Xios Plus	Orthophos DS
Xios XG	Pan Elite
Xios Scan	
Xios AE	
Orthophos XG (new plugin ¹)	
Orthophos XG 3D	
Orthophos SL, Orthophos S, Orthophos E	
Galileos	
FaceScan	

1. The OP XG plugin is available for download in the download section.

- Note: Individual product installation requirements and the document "Additional information on the operation with Sidexis 4" (download, 6524776)
- Sidexis 4 is only tested and released with the current PC-SW components.
Update all PC-SW components to the latest version (may require device software update and, if necessary, re-calibration)

Integration is possible for the cameras listed.

supported
SiroCam 3, C, digital, AF, AF+, F
SiroCam UAF Plus

- Appropriate plugins of the business division GE must be installed to be able to use these functions.
- The plugins and Siucom Plus required for the operation of the cameras are available in the operator area for the treatment centers.
- The individual product installation requirements must be met.

Integration of Dentsply Sirona treatment centers:

Dental Unit Plugins are available for Sidexis 4. When connected to Dentsply Sirona treatment centers, they enable

- Intraoral camera connection
- Additional display of the ApexLocator on a monitor*
- Automatic time synchronization*
- Patient-specific chair positions*
- Automatic head rest adjustment*

The plugins are available for download in the dealer area under the "Product Info -> Treatment Centers" menu item.

*not available for all Dentsply Sirona treatment centers

Other devices

- The Dürr IO scanner can be used via a new plugin from Dürr (new version of Bridge).
- Intraoral cameras from third-party suppliers can be connected via WDM (Windows Driver Model) to Sidexis 4 (Version 4.1.3 or later). To do so, install the video module from the Sidexis 4 Setup.
- As of CEREC Version 4.4, the CEREC database can be re-connected to Sidexis 4 via a new interface.
- SyncServices from REINKE SOFTWARE ENGINEERING (RSWE) and Sidexis 4 are fully compatible from Version 4.1.3 or later. A new version of SyncServices is required for this (to be obtained via RSWE).

Sidexis updates that are not yet supported

The following Sidexis XG updates are not yet supported by Sidexis 4 version 4.3:

- REPORTER
- TSE version

Should the customer continue to need these functions, then Sidexis XG is the more suitable product for these customers.

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Sirona Dental Systems GmbH